





U.S. Food and Drug Administration



CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

[FDA Home Page](#) | [CDRH Home Page](#) | [Search](#) | [CDRH A-Z Index](#) | [Contact CDRH](#)



[510\(k\)](#) | [Registration](#) | [Listing](#) | [Adverse Events](#) | [PMA](#) | [Classification](#) | [CLIA](#)  
[CFR Title 21](#) | [Advisory Committees](#) | [Assembler](#) | [Recalls](#) | [Guidance](#) | [Standards](#)

[New Search](#)

[Back To Search Results](#)

**510(k) Premarket Notification Database**

<b>Device Classification Name</b>	Mouthguard
<b>510(K) Number</b>	K053580
<b>Device Name</b>	DOCTOR'S NIGHTGUARD
<b>Applicant</b>	DENTAL CONCEPTS LLC. 555 Thirteenth Street, Nw Washington, DC 20004
<b>Contact</b>	Howard M Holstein
<b>Classification Product Code</b>	MQC
<b>Date Received</b>	12/22/2005
<b>Decision Date</b>	03/03/2006
<b>Decision</b>	Substantially Equivalent (SE)
<b>Classification Advisory Committee</b>	Dental
<b>Review Advisory Committee</b>	Dental
<b>Statement/Summary/Purged Status</b>	Summary Only
<b>Summary</b>	Summary
<b>Type</b>	Traditional
<b>Reviewed By Third Party</b>	No
<b>Expedited Review</b>	No

Database Updated 4/05/2007

[CDRH Home Page](#) | [CDRH A-Z Index](#) | [Contact CDRH](#) | [Accessibility](#) | [Disclaimer](#)  
[FDA Home Page](#) | [Search FDA Site](#) | [FDA A-Z Index](#) | [Contact FDA](#) | [HHS Home Page](#)

Center for Devices and Radiological Health / CDRH



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 26 2003

Dental Concepts LLC  
C/O Mr. Michael Lesser  
Medical Device Consultant, Incorporated  
49 Plain Street  
North Attleboro, Massachusetts 02760

Re: K024261

Trade/Device Name: Bite Plate  
Regulation Number: None  
Regulation Name: Dental Protector  
Regulatory Class: Unclassified  
Product Code: MQC  
Dated: March 5, 2003  
Received: March 6, 2003

Dear Mr. Lesser:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. Howard M. Holstein

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure.

**Indications for Use Statement**

510(k) Number (if known): K053580

Device Name: Doctor's® NightGuard™

**Indications for Use:**

The Doctor's NightGuard is indicated for protection against bruxism or nighttime teeth grinding. It is intended to reduce damage to the teeth and to prevent the noise associated with bruxing or grinding.

Prescription Use \_\_\_\_\_  
(Part 21 C.F.R. 801 Subpart D)

AND/OR

Over-The-Counter Use X  
(21 C.F.R. 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER  
PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan P. [Signature]

\_\_\_\_\_  
General Hospital  
Dental Devices

K053580

Page \_\_\_\_ of \_\_\_\_

K053580

MAR 3 2006

510(k) SUMMARY

**Dental Concepts The Doctor's® NightGuard™**

**Submitter's Name, Address, Telephone Number, Contact Person  
and Date Prepared**

Hogan & Hartson, LLP  
555 13<sup>th</sup> Street, N.W.  
Washington, D.C. 20004

Contact: Howard M. Holstein

Phone: (202) 637-5600

Facsimile: (202) 637-5910

Date Prepared: December 22, 2005

**Name of Device and Name/Address of Sponsor**

Doctor's® NightGuard™

Dental Concepts, LLC  
650 From Road  
Paramus, NJ 07652

Contact Person: Michael Lesser, President

Phone: (201) 225-2151

Facsimile: (201) 576-9780

**Common or Usual Name**

Dental Protector

**Classification Name**

Unclassified

### **Predicate Devices**

Dental Concepts BruxGuard  
Hollywood Products Mouth Peace  
GEM Scientific Products, Inc. Tension Reliever

### **Intended Use / Indications for Use**

The Doctor's NightGuard is indicated for protection against bruxism or nighttime teeth grinding. It is intended to reduce damage to the teeth and to prevent the noise associated with bruxing or grinding.

### **Technological Characteristics**

The Doctor's NightGuard is composed of a soft, formable clear upper material, made of ELVAX® resin, a copolymer of ethylene and vinyl acetate, and a hard occlusal base, which cushions the teeth. The base is composed of Elvaloy®, a copolymer of ethylene and methyl acrylate containing 9% methyl acrylate. When heated and then briefly cooled, the upper material can be molded to fit the user's upper teeth. The hard base prevents bite-through by users with moderate to severe nocturnal bruxing. The shock absorbing polymer material cushions the teeth on all sides.

### **Performance Data**

No performance data is required in support of this 510(k) notice.

### **Substantial Equivalence**

The Doctor's NightGuard is as safe and effective as Dental Concepts' BruxGuard. The two devices are physically the same. The labeling for the Doctor's NightGuard has been revised to make it suitable for OTC use. Because the two devices are the same, the Doctor's NightGuard possesses the same technological characteristics and principles of operation and a similar intended use as the BruxGuard predicate device. The Doctor's NightGuard also has similar intended uses and indications as the Hollywood Products Mouth Peace and the GEM Scientific Products, Inc. Tension Reliever, which were sold over the counter, like the NightGuard. Thus, the Doctor's NightGuard is substantially equivalent.